REF 42-9000

Denver Pleural Effusion Shunt

en	Pleural Effusion Shunt
fr	Shunt pour épanchement pleural
de	Pleuraeffusions-Shunt
it	Shunt per effusione pleurica
es	Shunt de efusión pleural
pt	Shunt de Efusão Pleural
nl	Pleurale effusie-shunt
sv	Pleuraexsudat-shunt
da	Pleuraekssudatshunt
fi	Pleuraeffuusiosuntti
no	Pleural effusjonsshunt
el	Αναστόμωση πλευρικής έγχυσης
ja	胸水シャント

USA Rx Only

STERILE EO

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© Copyright 2010, CareFusion Corporation,
or one of its subsidiaries. All rights reserved.

CareFusion
McGaw Park, IL 60085 USA

EC REP

CareFusion France 205 S.A.S.
F-44110 Châteaubriant

361-24009

• 2010-02

0123

For product inquiries or technical assistance:
pleurx-info@carefusion.com
www.carefusion.com/pleurx

Denver and Silique are trademarks and/or
registered trademarks of CareFusion
Corporation, or one of its subsidiaries.

CareFusion

Description

The Denver Pleural Effusion Shunt consists of a fenestrated pleural catheter, a flexible pump chamber containing two miter valves, and a fenestrated peritoneal catheter (**Figure 1**). The heart of the system is the valved pump chamber. This pump chamber is used by the patient or care giver to transfer the pleural fluid into the peritoneal cavity. The valves are made deliberately asymmetrical, causing them to be self-cleaning as the system is pumped. The second valve serves as a check that helps prevent reflux of liquid into the peritoneal catheter when the shunt is manually pumped.


Figure 1: Denver Pleural Effusion Shunt

The entire shunt is constructed from medical-grade silicone rubber. No connections at surgery are required; however, the peritoneal catheter may be shortened for optimal positioning in the peritoneal cavity. An integral barium sulfate stripe in the wall of both catheters permits visualization by x-ray or fluoroscopy to confirm proper positioning at the time of surgery.

Indications

The Denver Pleural Effusion shunt is indicated for use in patients with:

- chylothorax
- intractable aseptic pleural effusion

Caution: Benefits of instituting pleuro-peritoneal shunting in patients with malignant pleural effusion should be weighed against the potential hazards of dissemination of pleural fluid containing malignant cells.

Contraindications

The Denver Pleural Effusion Shunt is not designed for pleurovenous application.

Pleuro-peritoneal shunting is contraindicated in patients with infected pleural fluid and in patients suspected of having intra-abdominal infections such as diverticulitis.

Warnings

Never inject directly into the tubing portion of the shunt.

Do not use compressed air to open the valves or to dry the lumen of the shunt, since the valves or the pump body might balloon and burst.

Particulate matter may cause shunt occlusion, necessitating revision or replacement of the shunt. Resistance to compression of the pump chamber may be an indication that the peritoneal catheter or pump chamber is occluded. Should the pump chamber fail to refill, the pleural catheter may be occluded. (See Testing Shunt Function on page 5).

In general, shunting should be avoided if infection of any sort is present anywhere in the body. Infections of the shunt may require revision or removal of the shunt.

Careful placement of the catheters is important to avoid kinking, which may result in shunt malfunction.

Low-grade fever secondary to infusion of pleural fluid may occur and is self-limiting in the early postoperative period.

Cautions

For single use only. Re-use may result in a non-functional product or contribute to cross contamination.

Use sterile technique in all phases of the use of this product.

Benefits of instituting pleuro-peritoneal shunting in patients with malignant pleural effusion should be weighed against the potential hazards of dissemination of pleural fluid containing malignant cells.

Exercise extreme care to prevent the shunt from coming in contact with bare fingers, towels, drapes, talc, or any linty or granular surfaces. Silicone rubber is highly electrostatic and, as a result, attracts airborne particles and surface contaminants that could produce tissue reaction.

Use rubber-shod instruments when handling the shunt. Possible cuts or tears can occur if rubber-shod instruments are not used.

Exercise care when placing ligatures to avoid cutting or occluding the catheters. Do not use stainless steel ligatures.

Confirm shunt patency prior to insertion. Prime the shunt prior to insertion to eliminate large air bubbles from the pump chamber and catheters. (See Preparation of the Denver Pleural Effusion Shunt on page 2).

Potential complication of pleural effusion shunt implantation and use include, but may not be limited to the following: pneumothorax, laceration of lung or liver, wound infection, infection in the pleural space or peritoneal cavity, empyema, or subcutaneous tumor growth around the track of the shunt.

The patient or caregiver must be willing to pump the shunt on a regular basis (typically several minutes, three or four times a day). Without pumping, the shunt will not relieve symptoms of the pleural effusion.

The shunt can become occluded by proteins or other cellular debris. If occlusion occurs, the shunt must be partially or totally replaced to restore function.

Sterility

This product has been sterilized by ethylene oxide.

It is for single use only; do not resterilize.

CareFusion will not be responsible for any product that is resterilized, nor accept for credit or exchange any product that has been opened but not used.

Do not use if package is damaged.

To help prevent damage to sterile packaging, it is recommended that products be stored in the product box.

Preparation of the Denver Pleural Effusion Shunt

Caution: Confirm shunt patency prior to insertion. Prime the shunt prior to insertion to eliminate large air bubbles from the pump chamber and catheters.

Using sterile technique, immerse the entire shunt in sterile saline. Compress the pump chamber repeatedly until the system is filled with liquid and flow is established. **Note:** A small amount of crystalline sodium bicarbonate has been placed on the inside of the valve leaves to promote valve patency. The crystalline sodium bicarbonate will dissolve during preparation.

Expel any remaining air bubbles by elevating the peritoneal end of the shunt and alternately tapping and compressing the pump chamber while the entire pleural catheter remains immersed in the solution. Continue to pump the shunt body until all large air bubbles are removed.

Occasionally, the leaves of the silicone valve will adhere and the priming technique described above will fail to institute flow. This condition may be remedied by squeezing the miter valve edgewise while rolling the pump chamber at the site of the valve between the thumb and index finger (thus opening the valve), and then vigorously pumping liquid through the shunt. This maneuver assures patency of the shunt before insertion without damaging the valve mechanism.

Warning: Do not use compressed air to open the valves or to dry the lumen of the shunt, since the valves or the pump body might balloon and burst.

Suggested Operative Procedure

Caution: Use sterile technique in all phases of the use of this product.

Proper medical and surgical procedures are the responsibility of the surgeon. The appropriateness of any procedure must be based upon the needs of the patient.

This procedure may be performed under local anesthesia with preoperative sedation and analgesia, or general anesthesia. General anesthesia allows for positive pressure ventilation and may be helpful in limiting the risk of pneumothorax.

At the surgeon's discretion, the pleural catheter may alternatively be placed in the lower part of the pleural cavity instead of the upper part (**Figure 2**). If shorter catheter lengths are needed, the catheters may be cut to the desired length.


Figure 2: Placement

Place the patient in the supine position. Properly prep and drape the patient. Elevation of the patient's chest approximately 25° is useful but not essential.

Determine the proper intercostal space to insert the pleural catheter of the shunt. At the insertion site, make a transverse, inframammary incision (5-7 cm) in the mid-clavicular line. Carry the dissection down to the fascial layer (**Figure 3**).

Figure 3: Transverse Inframammary Incision

Use blunt dissection to create a subcutaneous pocket large enough to contain the pump chamber. The pocket should be inferior to the incision in such a position that the patient is able to see the bulge caused by the pump chamber (Figure 4).

Figure 4: Subcutaneous Pocket to Contain Pump

A 16 Fr. peel-away introducer will be used to place the pleural catheter of the shunt into the pleural cavity. Create a short subcutaneous tunnel, aimed superiorly and posteriorly from the lateral aspect of the incision (Figure 5).

Figure 5: Short Tunnel Created Above Rib Cage

The pleural catheter will be introduced into the pleural cavity at the most posterior and superior aspect of this tunnel. Pass the needle from the introducer set through the tunnel at the lateral aspect of the incision, just over the top of the rib (Figure 6).

Figure 6: Needle Placement

Use a syringe to aspirate pleural fluid. Insert the J-tip guidewire through the needle into the pleural space and remove the needle, leaving the guidewire in place (Figure 7).


Figure 7: J-Wire Insertion Through the Needle

Prepare the location for the insertion of the peritoneal catheter by making a 3 – 5 cm transverse, subcostal incision over the rectus muscle, on the same side as the pleural effusion, or just lateral to it. If the anterior rectus sheath is opened, split the rectus muscle, do not divide it. Expose the posterior sheath where the purse-string sutures are placed as shown (Figure 8). The alternative is to split the external and internal oblique muscles along the direction of their fibers. This exposes the transversalis fascia and peritoneum where the purse-string sutures are placed.


Figure 8: 3 to 5 cm Transverse Subcostal Incision

Pass the 16 Fr. peel-away introducer into the chest over the guidewire. Remove the dilator and guidewire, leaving the outer sheath in place. Pinch the sheath closed as the dilator and guidewire are removed to avoid inadvertent air entering the pleural space.

Pass a suction catheter through the sheath, removing as much pleural fluid as desired. Use of a suction catheter with an open arm allows the surgeon to control the applied negative pressure. Remove the suction catheter and insert the pleural catheter of the shunt through the sheath in an orientation parallel to the ribs and directed superiorly and posteriorly within the pleural cavity (Figure 9).

Figure 9: Insertion of Pleural Catheter of Shunt Through Sheath

Pinch the peel-away sheath closed during these maneuvers as necessary to avoid inadvertent air entering the pleural space. When the pleural catheter is in place, peel away the sheath and discard it.

Note: If a 16 Fr. peel-away introducer is not available, use a small clamp to tunnel superiorly and posteriorly as previously described. Then enter the pleural space, over the rib, with the clamp to establish access. Grasp the pleural catheter with the clamp which is used to carry the shunt catheter into the pleural space. Once in, divert the shunt posteriorly and superiorly unless the fluid is clearly located in a different anatomic orientation.

Pass the peritoneal catheter of the shunt subfascially from the thoracic to the abdominal incision, using a long clamp (Figure 10).

Warning: Careful placement of the catheters is important to avoid kinking, which may result in shunt malfunction.

Make a small incision inside the purse-string sutures and insert the catheter into the peritoneal cavity (Figure 10) under direct vision. Secure the purse-string sutures about the catheter taking care to achieve water tight closure without restricting the catheter lumen.


Figure 10: Subcutaneous Placement of the Shunt

Caution: Exercise care when placing ligatures to avoid cutting or occluding the catheters. Do not use stainless steel ligatures.

With the pump placed inferior to the chest incision, there is usually minimal discomfort associated with pumping. Secure the pump to the muscular fascia by placing several 3-0 monofilament sutures through the holes on the suture pad on either side of the pump chamber.

Compress the pump chamber several times to initiate fluid transfer. Irrigate both incision sites with antibiotic solution and close.

Postoperative Management

It is critical to ensure continuing shunt patency. Therefore, hourly compression of the pump chamber is recommended for the first 24 to 48 hours. This maneuver is best performed by the nursing staff for the first 12 to 24 hours following surgery. After that, the surgeon should be certain that the patient or caregiver is fully able to carry out instructions for pumping the shunt daily to maintain patency and transfer fluid.

Note: Each complete pump of the pump chamber will transfer approximately 2 ml of fluid. Steady pumping can result in the transfer of approximately 250 ml of fluid per 5 minute period of time. The number of pumps per day prescribed for the individual patient is determined by the surgeon, and is based on the surgeon's estimate of daily fluid accumulation.

Testing Shunt Function

Place the patient in the supine position. Grasp the sides of the pump chamber between the thumb and forefinger of one hand. Using a finger of the other hand, firmly compress the pump chamber. There should be little resistance to this compression. If there is resistance, then there may be an occlusion within the pump chamber or in the peritoneal catheter. If the pump chamber fails to refill readily after compression, the pleural catheter may be occluded.

Warning: Particulate matter may cause shunt occlusion, necessitating revision or replacement of the shunt. Resistance to compression of the pump chamber may be an indication that the peritoneal catheter or pump chamber is occluded. Should the pump chamber fail to refill, the pleural catheter may be occluded.

Replacement catheters and connectors are available for this purpose.

Ordering Information

REF	Description
42-9000	Pleural Effusion Shunt
10383-010	16 Fr. Introducer Set
42-3019	Nylon Connector (two per box) (for shunt revision)
42-2322	15.5 Fr. Fenestrated Catheter (27 cm)

Warranty

CareFusion warrants that this medical device is free from defects in both materials and workmanship. The above warranties are in lieu of all other warranties, either expressed or implied, including any warranty of merchantability or fitness for a particular purpose. Suitability for use of the medical device for any surgical procedure shall be determined by the user. CareFusion shall not be liable for incidental or consequential damages of any kind.

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