



Alaris® System Cleaning Products

The following cleaning products are recommended by CareFusion for use on the Alaris® System. Refer to the applicable Alaris® System Directions for Use for detailed cleaning instructions.

Cleaner	Dilution
10% bleach	1 part bleach to 9 parts water
70% isopropyl alcohol (IPA)	Not applicable
Compublend II (Base V with fragrance)	0.5 oz/gallon water
Aseptizyme	1 oz/gallon water
Clorox Wipes	Not applicable
Detergezyme	1 oz/gallon water
Dimension III	0.5 oz/gallon water
Dispatch	Not applicable
Expose II 256	0.5 oz/gallon water
Hibiclens	25.6 oz/gallon water
LpH Disinfectant Cleaner	0.5 oz/gallon water
Manu-Klenz	0.25 – 0.5 oz/gallon water
Maxima 128	1 oz/gallon water
Metrizyme	1 oz/gallon water
Mild detergent	per manufacturer's recommendation, as needed
Phenolic 256 DC	0.5 oz/gallon water
Professional Amphyl Hospital Bulk Disinfectant Cleaner	1.28 oz/gallon water
Staph-Attack	Not applicable
Super Sani-Cloth	Not applicable
Thymocide	Not applicable
Virkon	1.28 oz/gallon water
Warm water	Not applicable
Wexcide 128	1 oz/gallon water
Wexcide-Ready-To-Use	Not applicable