

Pyxis® Specimen Collection Verification

Identify and track specimens at the point of care



Designed to improve safety and efficiency during the laboratory specimen-collection process, Pyxis Specimen Collection Verification is an advanced wireless system that promotes positive patient identification and positive specimen identification for inpatients and outpatients at the point of care.

Pyxis Specimen Collection Verification requires barcode scanning of the patient wristband before blood is drawn or specimens are collected. The lightweight portable device includes a wireless Zebra Technologies printer to generate facility-specific barcode labels at the point of care. The system also can provide bedside access to patient demographics, allergy warnings, notifications, safety alerts, order verification and collection documentation. An optional collection status board provides a visual cue to enable Lab Managers to monitor collections and make adjustments as necessary, helping to improve and monitor productivity across the facility.

Pyxis®





Improve the specimen-collection process

- Track patient specimens from the point of care through the collection process
- Print facility-specific labels at point of care
- Implement safety checks to ensure right patient, order, test, container and time
- Document collection time, number of collections by personnel, incomplete draws and reasons for uncollected specimens automatically
- Support inpatient and outpatient collections drawn by either Nursing or Phlebotomy
- Optional collection status board enables Laboratory, Emergency Department, Respiratory and other Patient Care Managers, Physicians and Nurses to monitor the status of collections

Increase efficiency and help ensure patient safety through positive patient identification

- Verify any type of test order requiring collection of a specimen and LIS processing, including blood, urine, stool, sputum, etc
- Expedite stat draw by displaying prioritized test orders for collection personnel to view
- View patient-specific information and laboratory orders before visiting the patient
- Access incomplete collections audit trail to help track orders and assign or manage specific collections easily
- Scanning of patient and specimen information at the collection point helps ensure verification of positive patient identification throughout patient care

CareFusion
Reston, VA

carefusion.com

