

Imex Pocket-Dop II

Obstetric and Vascular Doppler



Obstetric and Vascular Versatility

Pocket-Dop II is a full-featured Doppler with interchangeable obstetric and vascular probes. Use the 2 or 3 MHz obstetric probes to monitor the fetal heartbeat. Use the 5 or 8 MHz vascular probes to screen patients at risk for stroke, peripheral arterial disease and other occlusive vascular conditions or to monitor difficult pulses and obtain systolic blood pressures quickly and easily.

Obstetric Examination

The Pocket-Dop II is a valuable aid in detecting pregnancy and monitoring the fetal heartbeat. The built-in speaker amplifies fetal sounds to reassure the expectant mother. A complete obstetric system should include a 3 MHz probe for early fetal heartbeat detection and a 2 MHz probe for use in later pregnancy or with larger patients. A 2 MHz waterproof probe is ideal for underwater labor or delivery.

Pocket-Dop II -
Obstetric
and Vascular
Screening and
Monitoring

Vascular Applications

Use the Pocket-Dop II to compare simple Doppler blood pressures in the arm and ankle to help diagnose arterial disease in the lower extremities. Or quickly screen for lower extremity vascular conditions such as valvular incompetence.

Emergencies and Challenging Clinical Situations

The Pocket-Dop II can be used to accurately monitor pulses and systolic blood pressures, even in noisy environments; confirm blood flow to wound sites, locate vessels for catheterization, withdrawal, injection or IV therapy; clearly hear low systolic blood pressures with patients who are diabetic, elderly or in shock.



System Configurations

Pocket-Dop II with:

2 MHz obstetric probe	P200
3 MHz obstetric probe	P300
5 MHz vascular probe	P500
8 MHz vascular probe	P800
Two 2 MHz obstetric probes	P220
2 & 3 MHz obstetric probes	P230
2 MHz obstetric & 5 MHz vascular probes ...	P250
2 MHz obstetric & 8 MHz vascular probes ...	P280
Two 3 MHz obstetric probes	P330
3 MHz obstetric & 5 MHz vascular probes ...	P350
3 MHz obstetric & 8 MHz vascular probes ...	P380
Two 5 MHz vascular probes	P550
5 & 8 MHz vascular probes	P580
Two 8 MHz vascular probes	P880
2 MHz waterproof obstetric probe	P2WP
No probe	PX00

Accessories

Probes:

2 MHz obstetric probe	T200
3 MHz obstetric probe	T300
5 MHz vascular probe	T500
8 MHz vascular probe	T800
2 MHz waterproof probe	WP20

Rollstands and Rollstand Accessories:

With the bracket	ST1
With bracket & storage basket	ST2
Storage basket for rollstand	BSKT
Locking bracket only	B100

Chargers:

110 volt	C610
220 volt	C619I
Battery	C620

Technical specifications:

Weight:	0.285 kg (10 oz.)
Dimension:	3.17 x 5.71 x 10.79 cm
.....	(1.25 x 2.25 x 4.25 in)
Power Source: ...	2 replaceable AA NiMH batteries
.....	110 volt or 220 volt battery charger
Warranty: 1 year labor and 5 years parts against	manufacture defects

This product is not CE marked, and not available for sale in the EU.

© 2011 CareFusion Corporation or one of its subsidiaries. All rights reserved
Lit. No. M80616 Rev 01 (2011/01)

CareFusion
Middleton, WI

carefusion.com

For more information, please contact your CareFusion
Distributor or call 877-842-7970, 303-431-9400 (Intl.)
or fax 303-431-0429

