

Imex Pocket-Dop 3

Obstetric and Vascular Doppler



Interchangeable, Convenient Probes

The Imex Pocket-Dop 3 is a valuable aid in detecting pregnancy and monitoring the fetal heartbeat. The built-in speaker amplifies fetal heart sounds to reassure the expectant mother. The Pocket-Dop 3 features interchangeable probes: the 3 MHz probe is for early fetal heartbeat detection and the 2 MHz probe for use in later pregnancy or with larger patients. The small "footprint" of these probes enables them to be easily angled down behind the pubic bone for optimum detection of early fetal heartbeats. When not in use, the probe is attached to the side of the unit by a unique magnetic holder. Obstetric exams are quick, easy and reliable with the Pocket-Dop 3.



Pocket-Dop 3 -
Obstetric Screening
and Monitoring

Accuracy of Autocorrelation

The Pocket-Dop 3 features autocorrelation which is a pattern recognition technique that quickly acquires and accurately tracks the complex signal of the fetal heart rate. This reduces the time wasted waiting for the system to acquire the best signal to read the fetal heart rate. The system averages the heart rate on command and then shows it in a high-contrast display panel. After three minutes of non-use, the unit automatically shuts off.

Smart Recharging System with Convenient Wall Mount

The Pocket-Dop 3 comes with a convenient recharging stand which can be placed on the counter or mounted to a wall. This smart recharging system won't overcharge the Doppler which ensures maximum battery life. The Pocket-Dop 3 is always handy and ready for use.



System Configurations

Pocket-Dop 3 with:

2 MHz obstetric probe	J200
Two 2 MHz obstetric probes	J220
3 MHz obstetric probe	J300
2 & 3 MHz obstetric probes	J230
Two 3 MHz obstetric probes	J330
No probe	JX00

Accessories

Probes:

2 MHz obstetric probe	H200
3 MHz obstetric probe	H300

Chargers:

110 volt	C640
220 volt	C640I
Charging Stand	C641
Battery	C622
Coiled Cord	CB0046

Technical Specifications

Dimensions:	12.7x8.38x2.5 cm (5x3.3x1 in.)
Weight:	0.285 kg (10 oz.)
Power Source:	
.....	8.4 volt NiMH rechargeable battery
.....	110 volt recharger
.....	220 volt recharger
Warranty: 5 years parts and 1 year labor against manufacture defects	

This product is not CE marked, and not available for sale in the EU.

© 2011 CareFusion Corporation or one of its subsidiaries. All rights reserved
Lit. No. M80881 Rev 01 (2011/01)

CareFusion
Middleton, WI

carefusion.com

For more information, please contact your CareFusion
Distributor or call 877-842-7970, 303-431-9400 (Intl.)
or fax 303-431-0429

