

Infant Flow[®] SiPAP System

Plus specifications



Modes
CPAP
CPAP with apnea/LBR
BiPhasic
BiPhasic with apnea/LBR

Alarms
High pressure
Low pressure
High FiO ₂
Low FiO ₂
Low battery
Low breath rate

Monitored parameters	
nCPAP	0 to 11 cmH ₂ O
MAP	0 to 11 cmH ₂ O
PIP	0 to 11 cmH ₂ O
PEEP	0 to 11 cmH ₂ O
%O ₂	21% to 100%

Warranty
The Infant Flow SiPAP System is warranted to be free from defects in material and workmanship and to meet published specifications for one year from date of shipment.


Controls	
nCPAP/pressure low flow meter	0 to 15 L/min
Pressure high flow meter	0 to 5 L/min
%O ₂	21% to 100%
Inspiratory time	0.1 to 3.0 seconds
Rate	1 to 54 bpm
Low breath rate interval	10 to 30 seconds (in 5 second intervals)
Manual breath	


Power supply
100 to 230 VAC; consumption 50 VA maximum
12 VDC (sealed lead acid battery); 2 hour life with full charge

Dimensions/weight (driver only)	
Dimensions	10.25" x 15" x 9.25" (26 cm x 38 cm x 23.5 cm)
Weight	19.5 lbs (8.8 kg)

Monitor screen	
Dimension	4.4" x 3.3" (11 cm x 8 cm)
Resolution	320 x 240 pixels, 1/4 VGA with auto-dimming backlight



 CareFusion
22745 Savi Ranch Parkway
Yorba Linda, CA 92887
800.231.2466 toll-free
714.283.2228 tel
714.283.8493 fax

 CareFusion Germany 234 GmbH
Leibnizstrasse 7
97204 Hoechberg
Germany
+49 931 4972-0 tel
+49 931 4972-423 fax



CareFusion
Yorba Linda, CA

carefusion.com

