



Switch-Blade[®] reusable laparoscopic scissors

Snowden-Pencer[®] laparoscopic products




Switch-Blade[®] reusable scissors with reusable handles and disposable scissor tips help provide superior cutting performance. The click-fit scissor tips make them easier to assemble by providing both audible and tactile feedback to confirm proper assembly.



Based on the next-generation Snowden-Pencer® ring handle design, the Switch-Blade® handle sets a new standard in comfort for laparoscopic instruments and accommodates both traditional precision-style grip or palm grip applications.

Assemble in three simple steps: thread, spread, click and then cut.




Thread
tip to shaft






Spread
the handle open



Click
to engage handle to tip



Catalog no.	Description	Qty
89-5103	Disposable curved METZ scissor tips 	10/box
89-5303	Disposable mini METZ scissor tips 	10/box
89-5115	Disposable straight scissor tips 	10/box
SP90-1053	Reusable handle and shaft, 36cm, insulated	1 each
SP90-1253	Reusable handle and shaft, 45cm, insulated	1 each



Click-fit tip attachment for easy assembly



Extended insulation helps prevent the exposure of any conductive metal surfaces

For more information or to place an order, please contact your V. Mueller® Snowden-Pencer® sales representative, call 800.323.9088 or visit carefusion.com/thinkblue.

© 2010 CareFusion Corporation or one of its subsidiaries. All rights reserved. Snowden-Pencer and Switch-Blade are registered trademarks of CareFusion Corporation or one of its subsidiaries. VM1073

CareFusion
Tucker, GA

carefusion.com

